Express Mail Label No.: ER833080649US

PCT/KR2003/000544 Docket No.: 59520-00009

10/521675 DT01 Rec'd PCT/PTS 18 JAN 2005

R. Comme

In the claims:

Please amend claim 6 and insert new claim 8 as shown below.

io Mad Stor 17

- 1. (Previously presented) Immunoglobin A polypeptide for diagnosing Diabetic retinopathy wherein said polypeptide is selected from the group consisting of a polypeptide sequence SEQ ID NO: 1 and a peptide fragment of the polypeptide sequence SEQ ID NO:1.
- 2. (Previously presented) The protein of claim 1 wherein the peptide fragment comprises peptide sequence SEQ ID NO: 2.
- 3-5. (Canceled).
- 6. (Currently amended) A method for diagnosing diabetic retinopathy, comprising:
 a) treating the antibody against polypeptide of claim 1 or 2 with a blood sample and an
 peroxidase, alkaline phosphatase or biotin conjugated-anti-Immunoglobulin A protein; and
 b) measuring optical density of the compound,
 wherein diabetic retinopathy is diagnosed when the measured value represents optical density
 (ELISA value) lower than normal one.
- 7. (Canceled).
- 8. (New) A method for diagnosing diabetic retinopathy, comprising:

A STATE OF THE STA

Express Mail Label No.: ER833080649US

PCT/KR2003/000544 Docket No.: 59520-00009

a) treating the antibody against polypeptide of claim 2 with a blood sample and an peroxidase, alkaline phosphatase or biotin conjugated-anti-Immunoglobulin A protein; and

b) measuring optical density of the compound,

wherein diabetic retinopathy is diagnosed when the measured value represents optical density (ELISA value) lower than normal one.